

# Determination in shellfish by liquid chromatography coupled to accurate mass spectrometry

Piva, E. et al. 2022 <https://doi.org/10.1002/dta.3282>



“...a rise in emerging chemicals that can’t be monitored with conventional methods, means we should explore alternative monitoring techniques.”

**Tarun Anumol, Ph.D.**

Global Director, Environmental Market, Agilent Technologies

## Study – Analyzed PFAS levels in mussel, clam and oyster samples



### Study



Used liquid chromatography and mass spectrometry (LC/MS) with Agilent 1290 Infinity II LC coupled to a 6546 quadrupole-time-of-flight mass spectrometer



LC/Q-TOF method used a library of 150 PFCA and PFSA compounds for non-targeted analysis



### Results



Every sample showed detectable PFAS



Of 12 PFAS detected, 7 PFAS were quantified (0.03-0.57 ng/g)



PFAS higher in Mediterranean than Pacific or Atlantic clams



N-MeFOSA, N-EtFOSA, N-MeFOSAA identified – possible new PFAS precursors



LC/MS/MS and LC/Q-TOF used for improved PFAS detection




### Conclusions



LC/Q-TOF method is effective in detecting linear and branched PFAS in shellfish at ng/g concentration



Method can also identify PFAS precursors e.g. 6:2 FTS and N-MeFOSAA in Mediterranean clams and mussels



Read the full Expert Insight, here



### Agilent Equipment Used

Agilent 6546 LC/Q-TOF and 1290 Infinity II HPLC



**WILEY** Analytical Science

# Perfluorooctane Sulfonate, Perfluorooctanoic Acid, and Hexafluoropropylene Oxide Dimer Acid (GenX) in a Benthic Fish

Hassell, K.L. *et al.* 2019 • <https://doi.org/10.1002/etc.4640>



“We are still at the tip of the iceberg in terms of information on PFAS occurrence and toxicity studies. Information is available for only a handful of them, but it’s prudent to monitor as many as possible now for baseline levels.”

**Tarun Anumol, Ph.D.**

Global Director, Environmental Market, Agilent Technologies



**Study on benthic fish on the Werribee River, Australia looked at:**  
PFAS depuration and accumulation • PFAS contaminations



## Study



40 male and 8 female adult blue spot gobies



Uniform food ingestion rates



11 week experiment



Ultra-high performance LC/MS/MS used to analyze body samples



Low limits of reporting for:

- PFOA
- Linear PFOS
- Linear+branched PFOS
- GenX



## Results



Accumulation of PFOA, linear PFOS and linear+branched PFOS



Steady-state concentration reached after 14 days



Depuration rate faster for PFOA compared to PFOS



Indication of different depuration kinetics for PFOS isomers



GenX (PFOS replacement) did not accumulate



## Research priorities



Specific PFOS isomers (branched and linear)




Tissue-specific distribution



Newer PFAS replacements



Organ-specific half-life data



Read the full Expert Insight, here



## Agilent Equipment Used

Agilent 6495 LC/TQ



**WILEY** Analytical Science