


Next-Gen CAR-T: High-Throughput Screening Accelerates Discovery


Revolutionizing antibody discovery and CAR construct optimization using automated high-throughput screening cytometry

The Problem and the High-Throughput Screening (HTS) Solution



CAR-T Therapy: Revolutionary but Bottlenecked


Engineered T cells targeting cancer-specific antigens



✓ FDA-approved for B-cell malignancies

✓ Durable remissions in patients with refractory malignancies


But...Traditional Development Is Slow



HTS: The Game Changer




Automated HTS Cytometry Transforms Every Stage



Traditional workflow

- Limited constructs per week
- Single parameter
- Days for results
- Manual processing

VS.



HTS Cytometry

- 100+ constructs/day
- 20+ parameters simultaneously
- Real-time kinetic data
- Automated analysis

Result



Vastly higher throughput + dramatically faster timelines



Significantly shortened development cycles

QC: Quality control; CAR-T: Chimeric antigen receptor T cell


Antibody Discovery and CAR Design Acceleration

Antibody Discovery: From Months to Weeks



The Challenge: Find optimal scFv with:

- ✓ High specificity
- ✓ Right affinity
- ✓ Minimal off-target effects
- ✓ Stable expression



Traditional Approach



Phage display



Individual clone testing




Low-throughput validation




Timeline: Many months


HTS Cytometry Approach



Parallel screening of antibody libraries in functional CAR-T format



Test hundreds of scFv candidates simultaneously



Measure binding + activation + function in one assay



Screen in relevant T cell context, not just binding



Timeline: Weeks


Key Advantage




Test scFvs directly as CAR constructs—see real functional impact immediately

CAR Design Optimization: Systematic Screening at Scale

CAR Components to Optimize:



Traditional: Test limited configurations (e.g., CD28 vs. 4-1BB)



HTS: Screen 40+ signaling domains in a parallel combinatorial library

Example Success: Novel Domain Discovery

Pooled CAR screening identified BAFF-R costimulatory domain



Enhanced cytotoxicity vs. standard 4-1BB



Better outcomes in xenograft models



Improved persistence under chronic antigen stimulation












From weeks of sequential testing to days of parallel profiling

scFv: Single-chain variable fragment









QC at Every Stage—Accelerated

Integrated QC Throughout Development









Use Case 1: Functional Validation (CD19-Targeted CAR-T)

Platform	Multiplexed Measurements (24–48 hours)			
<div> + </div> <div>Incucyte® + iQue® combined workflow</div>	<div></div> <div>Real-time tumor killing kinetics (Incucyte®)</div>	<div></div> <div>CAR-T activation (CD25 and CD69) (iQue®)</div>	<div></div> <div>Cytokine secretion (IFN-γ and TNF-α) (iQue®)</div>	<div></div> <div>T cell proliferation and phenotype (iQue®)</div>
Setup	Results			
Anti-CD19 CAR-T variants vs. CD19 + Ramos (target) and CD19-Jurkat (control) cells	CD: Cluster of differentiation; IFN-γ: Interferon-gamma; TNF-α: Tumor necrosis factor; EC50: Half-maximal effective concentration			
	<div></div> <div>Confirmed antigen-specific killing</div>	<div></div> <div>Quantified EC50 for each construct</div>	<div></div> <div>Minimal off-target activation</div>	
	vs.			
	Traditional: Multiple separate assays over extended timeframes			

Use Case 2: Exhaustion Profiling Under Chronic Stimulation

Challenge	iQue® Human T Cell Exhaustion Kit tracks:		
 Will CAR-T cells maintain function long-term?	 Exhaustion markers (PD-1, TIM-3, and LAG-3) expression	 Proliferative potential	 Cytokine production capacity (IFN-γ and TNF-α)
HTS Solution	Results		
 Continuous antigen challenge assay (10 days)	PD-1: Programmed cell death protein 1; TIM-3: T cell immunoglobulin and mucin-domain containing-3; LAG-3: Lymphocyte activation gene		
	 Identify CAR designs maintaining function under pressure	 Traditional: Single endpoint only	 HTS: Continuous monitoring = predictive data

Use Case 3: Safety Profiling—"On-Target, Off-Tumor" Risk


Critical Challenge	Test Panel			
 HER2 CAR-T toxicity in solid tumors	AU565 (high HER2)	MDA-MB-231 (low HER2)	MDA-MB-468 (negative)	+ HER2 CAR-T variants
Problem	Real-Time Incucyte® Monitoring			
 Attacks healthy tissue with low HER2 expression	 Strong killing of high HER2	 Minimal activity vs. low HER2	 No activity vs. negative	
HTS Solution	 Speed: Results in days vs. extended traditional timelines	 Scale: Screen 10+ CAR variants × 5+ cell lines in one 384-well experiment		
 3D spheroid specificity screen	Outcome	HER2: Human epidermal growth factor receptor 2		
	Identify tumor-selective, safer CAR designs before clinical testing			

The Integrated Platform and Impact

iQue® HTS Cytometry Platform



20+ parameter
flow cytometry



384-well
throughput



Cell phenotype + secreted
proteins (same well)



Automated ForeCyt®
software analysis



Incucyte® Live-Cell Analysis




Real-time kinetic
imaging (days-weeks)



96/384-well
automated acquisition



Label-free and
fluorescent detection



Continuous monitoring
in the incubator

Applications Across the CAR-T Pipeline



Discovery

- Antibody screening
- CAR construct libraries
- Specificity profiling



Optimization

- Killing kinetics
- Activation profiles
- Cytokine production



Expansion and QC

- Proliferation tracking
- Phenotype maintenance
- Exhaustion monitoring



Manufacturing


- Product potency
- Identity/purity
- Release testing

The Bottom Line: HTS Delivers



Speed

Dramatically faster
development timelines
(shortening months)



Scale

Vastly higher through-
put (100s of conditions
simultaneously)



Depth

20+ parameters vs.
1-3 traditional (richer
biological insight)



Predictive power

Real-time kinetics
(better clinical
translation)



Safety

Comprehensive
specificity profiling
(reduces clinical risk)

Result



More effective CAR-T therapies, reaching patients faster

WILEY

SARTORIUS

Further Resources



White paper:

[Phenotypic and Functional Characterization of CAR-T Cells with Advanced Flow Cytometry and Live-Cell Analysis](#)



Webinars:

• [CAR-T: Why Not Me?](#) • [The Importance of Immune Profiling in CAR-T Therapies](#)



Key Publications:

• [Wang et al. \(2021\) - High-throughput image cytometry for CAR-T. Cytometry Part A](#)

• [Sarikonda et al. \(2021\) - Best practices for CAR-T flow cytometry. Cytometry Part B](#)